


Test Case ID	Feature	Test Case Description	Expected Result	 Status	Device
TC01	Vitals Form	Open Record Vitals screen	Vitals form opens correctly	PASS	Mobile
TC02	Temperature	Enter temperature 36.5	Value accepted and displayed	PASS	Mobile
TC03	Blood Pressure	Enter BP 120/80	BP saved correctly	PASS	Mobile
TC04	Heart Rate	Enter heart rate 72	Value saved correctly	PASS	Mobile
TC05	Respiration Rate	Enter respiration rate 16	Value accepted	PASS	Mobile
TC06	SpO2	Enter oxygen level 98%	Value saved correctly	PASS	Mobile
TC07	Notes	Enter notes in notes field	Notes displayed correctly	PASS	Mobile
TC08	Save Vitals	Save vitals form	Data saved without error	PASS	Mobile
TC09	Review Vitals	Open saved vitals	Data appears correctly	PASS	Mobile
TC10	Edit Vitals	Edit vitals values	Updated values saved	PASS	Mobile
TC11	Timestamp	Check created timestamp	Correct time displayed	PASS	Mobile
TC12	Timestamp Update	Edit form and save	Updated timestamp appears	PASS	Mobile
TC13	Visit Note	Create Visit Note	Visit note saved successfully	PASS	Mobile
TC14	Visit Note Review	Open saved Visit Note	Notes appear correctly	PASS	Mobile
TC15	Visit Note Edit	Edit Visit Note	Updated notes visible	PASS	Mobile
TC16	AMPATH Form	Open AMPATH form	Form loads correctly	PASS	Mobile
TC17	AMPATH Save	Enter data and save	Form saved successfully	PASS	Mobile
TC18	AMPATH Edit	Edit AMPATH form	Updated values saved	PASS	Mobile
TC19	Date Picker	Tap date picker field	Calendar opens	PASS	Mobile
TC20	Date Selection	Select date from picker	Date appears in field	PASS	Mobile
TC21	Notes Formatting	Enter multi-line notes	Line breaks displayed	PASS	Mobile

Test Case ID	Feature	Test Case Description	Expected Result	 Status	Device
TC22	Empty Field	Save form with empty field	Validation error appears	FAIL	Mobile
TC23	Invalid Input	Enter text in numeric field	Error message shown	PASS	Mobile
TC24	UI Layout	Check form UI on mobile	Layout fits mobile screen	PASS	Mobile